



# *Product Catalogue*

EuniceMed Health Care



footcare / compression stocking / orthopedic support





# *Comfort Plus* *Corporation* EuniceMed Health Care

Comfort Plus Corporation, as an expert to design, develop, and manufacture the best quality of insoles and footcare, compression stocking, and orthopedic support products in Taiwan. Our research and development team has the professional technique capability, requisite knowledge, enthusiasm, and imagination to manage the constantly evolving technology and create new products to increase the market share. We have accumulated abundant experience and established stable and faithful long-term business relationship with our clients. Products are 100% made in Taiwan and promote to international global medical markets. All our products have been granted FDA & CE & ISO 13485 certificates and are selling well to global medical and healthcare markets.


EuniceMed of Comfort Plus Corporation continues to strive towards achieving its vision of leading in the design and manufacture service industry, developing the most advanced medical products one step at a time. We provide the best quality and reasonable price to look for distributors of the world.

We are also looking for our best OEM/ODM partnership and willing to be your best support to build up your own brand. With our professional technology and innovative spirit, we devoted ourselves to offer competitive price and high quality products to our clients. We believe a win-win situation especially for international business cooperation. To achieve this goal, we keep continuous product innovation and provide high quality products and services with competitive price and prompt delivery.

Our products are including following series

## Orthopedic support series

With different materials and design such as knitting, Neoprene, and foam and designs for knee, shin, ankle, elbow, wrist, finger, thumb, waist and neck to provide better protection.



## Compression stockings series

Various yarns provide different styles and functions with gradual decreased pressure design. Improves blood circulation and prevents varicose veins and edema legs.



## Insoles and Footcare series

Insoles which are made of medical grade silicone & gel materials, absorb shock, prevent uncomfortable friction and release pressure from feet.





# Orthopedic support series

EuniceMed's orthopedic support products which gather the innovative and professional information for excellent protection and omnibearing performance with palmary comfort.



**Breathable:** Breathable material brings comfort.



**Far infrared ray of Nano bamboo charcoal:** Nano bamboo charcoal with far infrared ray helps to improve circulation and keep health.



**Enhanced support:** Additional stay to enhanced support and prevent injuring when move with support.



**Silicone/Gel pad:** Increase the pressure and light massage on the joint.



**Strengthen strap provides comfortable compression:** Helps to stabilize the joint and prevents injuring.



**Patent product:** Innovative product with patent provides customers excellent protection and omnibearing performance with palmary comfort.

# Orthopedic support



CPO-1304

## Elbow support with silicone pads

### | Feature |

- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.
- Gel pad provides extra compression to stabilize the injured area.
- Provides suitable compression on the elbow joint.

S	15 - 20 cm	6 - 8"
M	20 - 25 cm	8 - 10"
L	25 - 30 cm	10 - 12"
XL	30 - 35 cm	12 - 14"

### | Material |

Nylon, Spandex, Gel pad

**QTY : 1 PC**

### | Application |

- Epicondylitis
- Tennis / Golfer's elbow



CPO-1305

## Bamboo elbow support

### | Feature |

- To maintain the structure of elbow joint and provide suitable compression on the elbow joint.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of elbow joint.
- Protect the weak elbow.

S	15 - 20 cm	6 - 8"
M	20 - 25 cm	8 - 10"
L	25 - 30 cm	10 - 12"
XL	30 - 35 cm	12 - 14"

### | Material |

Nylon, Spandex, Bamboo fiber.

**QTY : 1 PC**

### | Application |

- Weak elbow joint
- Elbow pain



CPO-1401

## Wrist support

### | Feature |

- To maintain the structure of wrist joint and provide suitable compression on the wrist joint.
- Protect the weak wrist.

S	14.5 - 16.5 cm	5 ¾ - 6 ½"
M	16.5 - 18.5 cm	6 ½ - 7 ¼"
L	18.5 - 21.5 cm	7 ¼ - 8 ½"
XL	21.5 - 23.5 cm	8 ½ - 9 ¼"

### | Material |

Nylon, Spandex, Bamboo fiber

**QTY : 1 PC**

### | Application |

- Weak wrist joint
- Wrist pain



CPO-1402

## Wrist sleeve with gel pad

### | Feature |

- Protection of weak or injured wrist joint.
- Soft gel pad release the pain of carpal tunnel syndrome.
- Decrease the probability of injured repetition.

Right/ Left		
S	13 - 17 cm	5 - 6¾"
L	17 - 25 cm	6¾ - 9¾"

### | Material |

Nylon, Spandex, Gel pad

**QTY : 1 PC**

### | Application |

- Carpal tunnel syndrome





CPO-1603

### Bamboo knee support

**| Feature |**

- To maintain the structure of knee joint and provide suitable compression on the knee joint.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of knee joint.
- Protect the weak Knee.

S	28.5 – 33.5cm	11 ¼ – 13 ¼"
M	33.5 – 39 cm	13¼ – 15 ½"
L	39 – 44.5 cm	15 ½– 17 ½"
XL	44.5 – 50 cm	17 ½ – 19 ¾"

**| Material |**

Nylon, Spandex, Bamboo fiber

**QTY : 1 PC**

**| Application |**

- Weak knee joint
- Knee pain



CPO-1604

### Ankle and shin sleeve

**| Feature |**

- Provide suitable compression on legs, and improve the performance of muscle strength.
- Provides firm support help to increase the stability of ankle joint.
- Open-toe design

S	28 – 31cm	11 – 12¼"
L	31 – 35 cm	12 ¼– 13 ¾"

**| Material |**

Nylon, Lycra Soft

**QTY : 1 PR**

**| Application |**

- Tired legs and ankles
- Weak shin



CPO-1605

### Shin sleeve

**| Feature |**

- Provide suitable compression on legs, and improve the performance of muscle strength.
- Helps relieve pain from an injured or weak shin.
- Easy on and off.

S	28 – 31cm	11 – 12 ¼"
L	31 – 35 cm	12 ¼– 13 ¾"

**| Material |**

Nylon, Lycra Soft

**QTY : 1 PR**

**| Application |**

- Weak shin
- Tired leg



CPO-1606

### Knee support with silicone pad

**| Feature |**

- Bilateral spiral stays help to support knee joint.
- Circular silicone pad for patella stabilization.
- To maintain the structure of knee joint and provide suitable compression on the knee joint.

S	33 - 36cm	13-14"
M	36 - 38cm	14-15"
L	38 - 41cm	15-16"
XL	41 - 46cm	16-18"
XXL	46 - 52cm	18-20½"

**| Material |**

Nylon, Spandex, Silicone pad, Steel stay

**QTY : 1 PC**

**| Application |**

- Osteoarthritis
- Patella instability
- Weak knee joint
- Knee pain



# Orthopedic support



CPO-1609

## Knee support with spiral stays

### | Feature |

- Bilateral spiral stays help to support knee joint.
- To maintain the structure of knee joint and provide suitable compression on the knee joint.
- Protect the weak Knee.

S	33 - 36cm	13-14"
M	36 - 38cm	14-15"
L	38 - 41cm	15-16"
XL	41 - 46cm	16-18"

### | Material |

Nylon, Spandex, Steel stay.

### QTY : 1 PC

### | Application |

- Osteoarthritis
- Weak knee joint
- Knee pain



CPO-1613

## Knee support with gel pad

### | Feature |

- Circular gel pad for patella stabilization.
- To maintain the structure of knee joint and provide suitable compression on the knee joint.
- Protect the weak Knee.

S	28.5 - 33.5cm	11 ¼ - 13 ¼"
M	33.5 - 39 cm	13 ¼ - 15 ½"
L	39 - 44.5 cm	15 ½- 17 ½"
XL	44.5 - 50 cm	17 ½ - 19 ¾"

### | Material |

Nylon, Spandex, Gel pad.

### QTY : 1 PC

### | Application |

- Patella instability
- Weak knee joint
- Knee pain



CPO-1703

## Bamboo ankle support

### | Feature |

- Protects the weak ankle joint, and relieves pain from injured.
- Provide suitable compression on the ankle joint.
- Far infrared ray of Nano bamboo charcoal improves body circulation for health care of ankle joint.

S	15.5 - 21cm	6 - 8 ¼"
M	21 - 26 cm	8 ¼ - 10 ¼"
L	26 - 31 cm	10 ¼ - 12 ¼"
XL	31 - 36 cm	12 ¼ - 14 ¼"

### | Material |

Nylon, Spandex, Bamboo fiber.

### QTY : 1 PC

### | Application |

- Ankle pain
- Weak ankle joint



CPO-1704

## Ankle support with silicone pads

### | Feature |

- Bilateral silicone pads provide soft compression and stabilize the ankle.
- Provide suitable compression on the ankle joint.

S	15-20cm	6-8"
M	20-25cm	8-10"
L	25-30cm	10-12"
XL	30-35cm	12-14"

### | Material |

Nylon, Spandex, Silicone pad.

### QTY : 1 PC

### | Application |

- Ankle instability
- Ankle pain or weakness
- Ankle sprain





CPO-2301

### Tennis elbow support

**| Feature |**

- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.
- Silicone pad provides extra compression to stabilize the injured area.

One Size

**| Material |**

Neoprene, Nylon, Silicone pad.

**QTY : 1 PC**

**| Application |**

- Tennis elbow
- Golfer's elbow



CPO-2302

### Tennis/Golf elbow support

**| Feature |**

- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.
- Silicone pad provides extra compression to stabilize the injured area.
- Light weight design helps to increase the performance for sport.
- Easy to wear and adjustable.

One Size

**| Material |**

Neoprene, Nylon, Silicone pad.

**QTY : 1 PC**

**| Application |**

- Tennis elbow
- Golfer's elbow



CPO-2401

### Wrist and thumb support

**| Feature |**

- Provide mild compression and mild support to palms, thumb and wrist joint.
- Elastic hook and loop strap may be used to limit wrist joint.
- Provide basic support, relieve stiffness of hand or pain of wrist.

S	14.5 – 16.5cm	5 ¾ – 6 ½"
M	16.5 – 18.5 cm	6 ½ – 7 ¼"
L	18.5 – 21.5 cm	7 ¼ – 8 ½"
XL	21.5 – 23.5 cm	8 ½ – 9 ¼"

**| Material |**

Neoprene, Nylon, Spandex.

**QTY : 1 PC**

**| Application |**

- Stiffness hand
- Wrist weakness
- Wrist sprain



CPO-2402

### Thumb splint

**| Feature |**

- For De Quervain's syndrome use.
- Aluminum stay limited the range of motion of the thumb, and decrease painful and the probability of injured repetition.

S/M	14.5-18.5cm	5 ¾ – 7 ¼"
L/XL	18.5-23.5cm	7 ¼ – 9 ¼"

**| Material |**

Neoprene, Nylon, Spandex, Aluminum stay.

**QTY : 1 PC**

**| Application |**

- De Quervain's syndrome
- Thumb sprain



# Orthopedic support



## Enhanced wrist splint

### | Feature |

- For carpal tunnel syndrome use.
- Two aluminum stays enhance the stability of wrist, and decrease painful and the probability of injured repetition.

### Right/ Left

S	14.5 – 16.5cm	5 ¾ – 6 ½"
M	16.5 – 18.5 cm	6 ½ – 7 ¼"
L	18.5 – 21.5 cm	7 ¼ – 8 ½"
XL	21.5 – 23.5 cm	8 ½ – 9 ¼"

### | Material |

Neoprene, Nylon, Aluminum stay.

### QTY : 1 PC

### | Application |

- Carpal tunnel syndrome
- Wrist sprain
- Colles fracture



## Wrist wrap

### | Feature |

- Provide stability for wrist joint.
- Easy to wear and adjustable.

One Size

### | Material |

TPE, Nylon, Polyester.

### QTY : 1 PC

### | Application |

- Wrist pain
- Wrist weakness



## Wrist wrap

### | Feature |

- Provide stability for wrist joint.
- Easy to wear and adjustable.
- Breathable elastic material decrease muggy feeling.

One Size

### | Material |

Neoprene, Nylon, Spandex.

### QTY : 1 PC

### | Application |

- Wrist pain
- Wrist weakness



## Knee support

### | Feature |

- Easy to wear and adjustable.
- For knee joint stabilization.
- Comfortable for weak muscles and knee.

One Size

### | Material |

Neoprene, Nylon.

### QTY : 1 PC

### | Application |

- Knee weakness
- Knee pain





### Hinged knee brace

**| Feature |**

- Bilateral hinges help to stabilize knee joint.
- Opening design in the front to stabilize patella.
- Hook and loop strips provide comfortable compression and perfectly fit.

S	33 – 37cm	13 – 14 ½"
M	37 – 41 cm	14 ½ – 16"
L	41 – 45 cm	16 – 17 ¾"
XL	45 – 49 cm	17 ¾ – 19 ¼"

**| Material |**

Neoprene, Nylon, Spandex, Steel stay.

**QTY : 1 PC**

**| Application |**

- Weak or overstressed knee
- Osteoarthritis
- Mild strain / sprain of knee



### Knee brace with circular strip

**| Feature |**

- Bilateral spiral stays help to support knee joint.
- Opening design in the front with silicone pads provides soft compression to stabilize patella.
- Extra strip helps to stabilize the knee joint.

S	33 – 37cm	13 – 14 ½"
M	37 – 41 cm	14 ½ – 16"
L	41 – 45 cm	16 – 17 ¾"
XL	45 – 49 cm	17 ¾ – 19 ¼"

**| Material |**

Neoprene, Nylon, Silicone pad, Steel stay.

**QTY : 1 PC**

**| Application |**

- Weak or overstressed knee
- Osteoarthritis
- Mild strain / sprain of knee



### Knee support with open knee

**| Feature |**

- Easy to wear and adjustable.
- Opening design with silicone pad provides softness compression to stabilize patella.

One Size

**| Material |**

Neoprene, Nylon, Silicone pad.

**QTY : 1 PC**

**| Application |**

- Patella instability
- Knee weakness
- Knee pain



### Patella strap

**| Feature |**

- Easy to wear and adjustable.
- EVA pad provides extra compression on patellar tendon.
- Relieves pain from Chondromalacia Patella.

One Size

**| Material |**

Neoprene, Nylon, Spandex, EVA.

**QTY : 1 PC**

**| Application |**

- Chondromalacia Patella
- Patella instability



# Orthopedic support



CPO-2607

## Knee support with patella wrap

### | Feature |

- Silicone pad provides extra compression on patellar tendon.
- Spiral stays help to support knee joint.
- Breathable Neoprene provides comfortable wear.

S	33 – 37cm	13 – 14 ½"
M	37 – 41 cm	14 ½ – 16"
L	41 – 45 cm	16 – 17 ¾"
XL	45 – 49 cm	17 ¾ – 19 ¼"

### | Material |

Neoprene, Nylon, Steel stay, silicone pad.

### QTY : 1 PC

### | Application |

- Patella instability
- Knee weakness
- Knee pain



CPO-2608

## Knee support with open patella

### | Feature |

- Easy to wear and adjustable.
- Circular pad with open knee design for patella stabilization.
- Spiral stays help to support knee joint.

One Size

### | Material |

Neoprene, Nylon, Polyester, Steel stay.

### QTY : 1 PC

### | Application |

- Patella instability
- Knee weakness
- Knee pain



CPO-2701

## Ankle support

### | Feature |

- Easy to wear and adjustable.
- For ankle joint stabilization.
- Assist the rehabilitation of ankle sprain.

S	15.5 – 21cm	6 – 8 ¼"
M	21 – 26 cm	8 ¼ – 10 ¼"
L	26 – 31 cm	10 ¼ – 12 ¼"
XL	31 – 36 cm	12 ¼ – 14 ¼"

### | Material |

Neoprene, Nylon.

### QTY : 1 PC

### | Application |

- Ankle instability
- Ankle pain or weakness
- Ankle sprain / strain



CPO-2702

## Ankle support

### | Feature |

- Easy to wear and adjustable.
- For ankle joint stabilization.
- Helps relieve aching muscles and injured ankle.

S	15.5 – 21cm	6 – 8 ¼"
M	21 – 26 cm	8 ¼ – 10 ¼"
L	26 – 31 cm	10 ¼ – 12 ¼"
XL	31 – 36 cm	12 ¼ – 14 ¼"

### | Material |

Neoprene, Nylon

### QTY : 1 PC

### | Application |

- Ankle instability
- Ankle pain or weakness
- Ankle sprain / strain



CPO-2705

### Ankle support with Achilles protection pad

**| Feature |**

- Easy to wear and adjustable.
- Gel pad provides protection for Achilles tendon.
- For ankle joint stabilization.

**| Material |**

Neoprene, Nylon, Polyester, Gel pad.

**QTY : 1 PC**

**| Application |**

- Achilles tendinitis
- Ankle instability
- Ankle sprain / strain

S	20 – 26 cm	8 – 10 ¼"
L	26 – 32 cm	10 ¼ – 12 ½"



CPO-2706

### Ankle support

**| Feature |**

- Easy to wear and adjustable.
- Helps relieve aching muscles and injured ankle.
- For ankle joint stabilization.

**| Material |**

Neoprene, Nylon

**QTY : 1 PC**

**| Application |**

- Ankle pain or weakness
- Ankle instability
- Ankle sprain / strain

S	15.5 – 21cm	6 – 8 ¼"
M	21 – 26 cm	8 ¼ – 10 ¼"
L	26 – 31 cm	10 ¼ – 12 ¼"
XL	31 – 36 cm	12 ¼ – 14 ¼"



CPO-3101

### Foam cervical collar

**| Feature |**

- To stabilize the cervical spine in neutral position.
- velcro tape for easy adjustment.

**| Material |**

cotton yarn with F/I, PU

**QTY : 1 PC**

**| Application |**

- Cervical soreness
- Cervical sprain / strain

S	25.5 – 33cm	10 – 13"
M	33 – 40.5 cm	13 – 16"
L	40.5 – 48 cm	16 – 19"



CPO-4101

### Philadelphia collar

**| Feature |**

- Provides the stability and limited the range of motion of cervical vertebra.
- Anterior opening for trachea.
- Hook and loop closure for easy adjustment.

**| Material |**

Foam, PP, Nylon.

**QTY : 1 PC**

**| Application |**

- Cervical spondylosis
- Post-operative use for cervical care
- Cervical fracture

S	25.5 – 33cm	10 – 13"
M	33 – 40.5 cm	13 – 16"
L	40.5 – 48 cm	16 – 19"



# Orthopedic support



## Curved finger splint

### | Feature |

- For stabilizes of the dislocated finger joints or fracture of the phalanges .

S	7.5 cm
M	10.3 cm
L	14 cm

### | Material |

Aluminum, PE.

**QTY : 1 PC**

### | Application |

- Dislocation injury of phalanges
- Fracture of the phalanges



## Fold-over finger splint

### | Feature |

- For protecting injuries to the finger tips.
- For total finger immobilization.

S	4.5 cm
M	7 cm
L	8.5 cm

### | Material |

Aluminum, PE.

**QTY : 1 PC**

### | Application |

- Hyperextension and dislocation injuries of phalanges
- Collateral ligament injury of distal interphalangeal joint



## Fold-over finger splint

### | Feature |

- For protecting injuries to the finger tips.
- For total finger immobilization.

S	4.5 cm
M	7 cm
L	8.5 cm

### | Material |

Aluminum, PE.

**QTY : 1 PC**

### | Application |

- Hyperextension and dislocation injuries of phalanges
- Collateral ligament injury of distal interphalangeal joint
- Fracture of the phalanges



## Four prong finger splint

### | Feature |

- Protects injuries to finger tips.
- For total finger immobilization.

S	5 cm
M	6 cm
L	7 cm

### | Material |

Aluminum, PE.

**QTY : 1 PC**

### | Application |

- Hyperextension and dislocation injuries of phalanges
- Collateral ligament injury of distal interphalangeal joint
- Fracture of the phalanges





### Frog finger splint

**| Feature |**

- To hold and protect the injured finger joint in a desired alignment.

S	5.5 cm
M	7.5 cm
L	8 cm

**| Material |**

Aluminum, PE.

**QTY : 1 PC**

**| Application |**

- Dislocation injury of phalanges
- Fracture of the phalanges



### Baseball finger splint

**| Feature |**

- For stabilizes the fractures of the distal phalanges.
- Can be adjusted to amount of extension or hyperextension.

S	7 cm
M	10 cm
L	10.5 cm

**| Material |**

Aluminum, PE.

**QTY : 1 PC**

**| Application |**

- Dislocation injury of phalanges
- Fracture of the phalanges



### Plastic finger splint

**| Feature |**

- For stabilizes of the dislocated finger joints or fracture of the phalanges.
- To protect the injured finger.

	height	width	length
1	14cm	17cm	50cm
2	15cm	17cm	52cm
3	18cm	19cm	56cm
4	19cm	19cm	61cm
5	19cm	20cm	56cm
5.5	22cm	24cm	60cm
6	23cm	26cm	61cm
7	25cm	27cm	62cm

**| Material |**

ABS

**QTY : 1 PC**

**| Application |**

- Hyperextension and dislocation injuries of phalanges
- Collateral ligament injury of distal interphalangeal joint
- Fracture of the phalanges



### Air/gel ankle brace

**| Feature |**

- Half air and half gel design to bind ankle joint.
- The splints in both sides to stabilize injured ankle joint from inversion and eversion.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

**| Material |**

PP, PVC, Nylon, Gel.

**QTY : 1 PC**

**| Application |**

- Ankle sprain / strain
- Acute injuries of ankle joint
- Mild swollen ankle



# Orthopedic support



CPO-4702

## Adjustable air ankle brace

### | Feature |

- Adjustable air pad design provide perfectly fit on ankle joint.
- The splints in both sides to stabilize injured ankle joint from inversion and eversion.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

### | Material |

PP, PVC, Nylon, PV.

**QTY : 1 PC**

### | Application |

- Ankle sprain / strain
- Acute injuries of ankle joint
- Mild swollen ankle



CPO-4703

## Gel ankle brace

### | Feature |

- Cold compression therapy to ankle joint by iced insert pad.
- Insert pad could be took off and put in the refrigerator before use.
- The splints in both sides to stabilize injured ankle joint from inversion and eversion.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½

### | Material |

PP, Nylon, Gel, PV.

**QTY : 1 PC**

### | Application |

- Ankle sprain / strain
- Acute injuries of ankle joint
- Mild swollen ankle



CPO-6403

## Wrist support

### | Feature |

- Protection of weak or injured wrist joint.
- Provide suitable compression and stability for wrist joint.

One Size

### | Material |

Nylon, Rubber, Polyester

**QTY : 1 PC**

### | Application |

- Wrist pain
- Wrist weakness



CPO-6404

## Thumb splint

### | Feature |

- For De Quervain's syndrome use.
- Aluminum stay limited the range of motion of the thumb, and decrease painful and the probability of injured repetition.
- Knitting material decrease muggy feeling and more fitting to hand shape.

S	13 - 17 cm	5 - 6¾"
L	17 - 25 cm	6¾ - 9¾"

### | Material |

Cotton, Rayon, Rubber, Aluminum

**QTY : 1 PC**

### | Application |

- De Quervain's syndrome
- Thumb sprain





CPO-6202

### Lumbar support

**| Feature |**

- 9 inches height.
- 4 steel stays provide stable support of back.
- Outside strengthen strap increases perfectly fit and compression of the belt.

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"
XXL	109 - 119 cm	43 - 47"

**| Material |**

Polyester, Cotton, Rubber, Nylon, Steel stay.

**QTY : 1 PC**

**| Application |**

- Low back pain
- Lumbar spondylosis
- Spondylolisthesis
- Back strain



CPO-6204

### Breathable lumbar support

**| Feature |**

- 8 inches height.
- 4 steel stays provide stable support of back.
- Slim breathable mesh decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"

**| Material |**

Nylon, Polyester, Cotton, Rubber, Steel stay.

**QTY : 1 PC**

**| Application |**

- Low back pain
- Lumbar spondylosis
- Back strain



CPO-6206

### Magnetic back belt

**| Feature |**

- Magnet design in the waist and belly help to increase circulate.
- 6 inches height, well binding.
- Suitable for health care use.

S	58 - 68 cm	23 - 27"
M	68 - 80 cm	27 - 31 1/2"
L	80 - 92 cm	31 1/2 - 36"
XL	92 - 105 cm	36 - 41 1/4"

**| Material |**

Polyester, Nylon, Cotton, Rubber, Magnet.

**QTY : 1 PC**

**| Application |**

- Low back pain
- Back strain



CPO-6207

### Lumbar belt

**| Feature |**

- 6 inches height, perfectly fit.
- 2 plastic stays for health care use.
- Increase the abdominal pressure to promote the muscle strength in waist.
- Elastic strap provides maximum comfort, lightweight and

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"

**| Material |**

Cotton, Nylon, Rubber.

**QTY : 1 PC**

**| Application |**

- Low back pain
- Back strain



# Orthopedic support



CPO-6208

## 6" high lumbar support

### | Feature |

- 6 inches height, perfectly fit.
- 4 plastic stays provide slightly support in waist.
- Outside strengthen strap increases perfectly fit of the belt.

### | Material |

Cotton, Nylon, Polyester, Rubber, PE.

### QTY : 1 PC

### | Application |

- Low back pain
- Lumbar spondylosis
- Back strain

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"



CPO-6211

## Breathable lumbar support

### | Feature |

- 9 inches height.
- 4 steel stays provide stable support of back.
- Slim breathable mesh decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

### | Material |


Nylon, Polyester, Rubber, Steel stay

### QTY : 1 PC

### | Application |

- Low back pain
- Lumbar spondylosis
- Back strain

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"



CPO-6212

## Industrial back support

### | Feature |

- 8 inches height.
- 4 plastic stays provide stable support of back.
- Adjustable shoulder suspension straps.
- Outside strengthen strap increases perfectly fit and compression of the belt.

### | Material |

Nylon, Polyester, Rubber, Steel stay

### QTY : 1 PC

### | Application |

- Low back pain
- Back strain

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"





CPO-6213

### Breathable lumbar support

**| Feature |**

- 8 inches height.
- 2 plastic stays provide stable support of back.
- Slim breathable mesh decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

**| Material |**

Nylon, Polyester, Rubber, Cotton, PP

**QTY : 1 PC**

**| Application |**

- Low back pain
- Back strain

S	58 - 68 cm	23 - 27"
M	68 - 80 cm	27 - 31 1/2"
L	80 - 92 cm	31 1/2 - 36"
XL	92 - 105 cm	36 - 41 1/4"



CPO-6214

### lumbar support

**| Feature |**

- 9 inches height.
- 4 steel stays provide stable support of back.
- Outside strengthen strap increases perfectly fit and compression of the belt.

**| Material |**

Nylon, Polyester, Rubber, Steel stay

**QTY : 1 PC**

**| Application |**

- Low back pain
- Lumbar spondylosis
- Back strain
- Spondylolisthesis

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"



CPO-6223

### Breathable lumbar support

**| Feature |**

- 9 inches height.
- 4 plastic stays provide slightly support in waist.
- Slim breathable mesh decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

**| Material |**

Nylon, Polyester, PP

**QTY : 1 PC**

**| Application |**

- Low back pain
- Lumbar spondylosis
- Back strain

S	58 - 70 cm	23 - 27 1/2"
M	70 - 85 cm	27 1/2 - 33 1/2"
L	85 - 100 cm	33 1/2 - 39 1/2"
XL	100 - 110 cm	39 1/2 - 43 1/4"
XXL	109 - 119 cm	43 - 47"





CPO-6601

## Elastic knee support

### | Feature |

- Bilateral spiral stays provide stable support for knee joint.
- Breathable elastic material decrease muggy feeling.
- To maintain the structure of knee joint and provide suitable compression on the knee joint.

### | Material |

Polyester, Rayon, Steel stay

### QTY : 1 PC

### | Application |

- Patella instability
- Knee weakness
- Knee pain

<b>S</b>	30 - 34 cm	11 7/8 - 13 3/8"
<b>M</b>	34 - 38 cm	13 3/8 - 15"
<b>L</b>	38 - 42 cm	15 - 16 1/2"
<b>XL</b>	42 - 46 cm	16 1/2 - 18 1/8"
<b>XXL</b>	46 - 50 cm	18 1/8 - 19 3/4"



CPO-7202

## Lumbar support

### | Feature |

- 9 inches height.
- 4 steel stays provide stable support of back.
- Breathable mesh fabric decrease muggy feeling.
- Outside strengthen strap increases perfectly fit and compression of the belt.

### | Material |

Nylon, Polyester, PP, Steel stay.

### QTY : 1 PC

### | Application |

- Low back pain
- Lumbar spondylosis
- Back strain

<b>S</b>	58 - 70 cm	23 - 27 1/2"
<b>M</b>	70 - 85 cm	27 1/2 - 33 1/2"
<b>L</b>	85 - 100 cm	33 1/2 - 39 1/2"
<b>XL</b>	100 - 110 cm	39 1/2 - 43 1/4"



Patent  
CPO-7304

## Tennis/Golf elbow support

### | Feature |

- To tighten and support the injured muscles and tendons in tennis elbow and golfer's elbow injuries.
- Silicone pad provides extra compression to stabilize the injured area.
- Light weight design helps to increase the performance for sport.
- Easy to wear and adjustable.

### | Material |

Neoprene, Nylon, Silicone pad.

### QTY : 1 PC

### | Application |

- Tennis elbow
- Golfer's elbow

<b>S</b>	20 - 28 cm	8 - 11"
<b>L</b>	28 - 34 cm	11 - 13 1/2"





CPO-7401

## Wrist brace

### | Feature |

- For carpal tunnel syndrome use.
- Aluminum stay enhance the stability of wrist joint.
- Easy to wear and adjust.

### | Material |

Polyester, Nylon, Aluminum stay.

**QTY : 1 PC**

### | Application |

- Carpal tunnel syndrome
- Wrist sprain
- Colles fracture

<b>S</b>	14.5 – 16.5cm	5 3/4 – 6 1/2"
<b>M</b>	16.5 – 18.5 cm	6 1/2 – 7 1/4"
<b>L</b>	18.5 – 21.5 cm	7 1/4 – 8 1/2"
<b>XL</b>	21.5 – 23.5 cm	8 1/2 – 9 1/4"



CPO-7701

## Ankle brace

### | Feature |

- Assist the rehabilitation of ankle sprain
- Two plastic pads enhance the stability.
- For ankle joint stabilization.

### | Material |

Nylon, Polyester, PE

**QTY : 1 PC**

### | Application |

- Ankle instability
- Ankle sprain / strain

<b>S</b>	15.5 - 21 cm	6 - 8 1/4"
<b>M</b>	21 - 26 cm	8 1/4 - 10 1/4"
<b>L</b>	26 - 31 cm	10 1/4 - 12 1/4"
<b>XL</b>	31 - 36 cm	12 1/4 - 14 1/4"



Patent  
CPO-7704

## Toe splint for Hallux valgus

### | Feature |

- Bilateral light aluminum hinges serve as the prop of toes and keep them fixed.
- Detachable hinge on little toes can be applied if need.
- Relieves pain caused by hallux valgus and Reduces friction between toes.
- Breathable mesh fabric for comfort feeling.
- Can be used while walking.


### | Material |

Neoprene, Nylon, Aluminum, Polyester

**QTY : 1 PC**

### | Application |

- Hallux valgus.
- Tailor's bunion.
- Varus toe.



<b>S</b>	20 - 24 cm	8 - 9 1/2"
<b>L</b>	24 - 29 cm	9 1/2 - 11 1/2"





# Compression stockings series

*Medical compression stockings help fight blood pressure within the veins, mitigate or avoid aggravated of damage. Graduated compression design helps promote better blood flow and circulation, it can help the blood return and increase the function of circulation metabolic.*

## Features of Compression Stockings

- Design for German regulatory compression standard, ccl.1:18-21mmHg, ccl.2:23-32mmHg, and ccl.3:34-43mmHg. Graduated compression design helps promote better blood flow and circulation, reduce the edema and aching of lower extremity.
- Steady state of graduated compression helps legs in good status, effectively eliminate the edema and shape the legs.
- Lycra yarn has better elasticity and durability, it can maintain the efficacy of product and not be friable easily. The soft hand feel provides users for long time wearing comforts.
- Non-toxic yarn reduce the allergies, the yarn with the certificate of Oeko-Tex Standard 100 make consumers more reassurance.
- All of long-term standing people, teacher, salesperson, waitress, flight attendant, and person who has edema or tired, aching for lower limbs are suitable for ccl.1 compression stockings. Someone who is slight varicosity can use ccl.1 to avoid varicosity get worse. Moderate varicosity or the patient who was post surgical can use ccl.3 for recovery, and use ccl.2 as the long term maintenance.

## Applications of each compression level

### ccl.1:18-21mmHg

- Minor varicosity.
- Tired, aching legs.
- Minor ankle, leg and foot swelling.
- Reduce leg swelling and prevent varicosity during pregnancy.

### ccl.2:23-32mmHg

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

### ccl.3:34-46mmHg

- Severe varicosity.
- Severe edema.
- Post surgical use.



Graduated compression design helps promote better blood flow and circulation.



Integral knitting for heel and toe design provides comfortable and easy wearing.



Patent Lycra material to provide more soft touch and long time wearing.



Elastic band prevents slipping and binding.



Open toe design helps to check blood circulation.



Soft band with silicone dots prevents slipping and binding.



Special design with less elasticity for hip helps to improve the shape of the hip and the circulation of legs.

# Compression stockings

Compression level

## Class 1 (18-21 mmHg)



### Knee high compression stockings (Closed toe)

**| Application |**

- Minor varicosity.
- Tired, Aching Legs.
- Minor ankle, leg and foot swelling.
- Reduce leg swelling during pregnancy.

**| Material |**

Lycra, Nylon.

Size 1 – 5 (Refer to P.29)

Color  



### Mid thigh compression stockings (Closed toe)

**| Application |**

- Minor varicosity.
- Tired, Aching Legs.
- Minor ankle, leg and foot swelling.
- Reduce leg swelling during pregnancy.

**| Material |**

Lycra, Nylon.

Size 1 – 5 (Refer to P.29)

Color  



### Thigh high compression stockings with silicone lace band (Closed toe)

**| Application |**

- Minor varicosity.
- Tired, Aching Legs.
- Minor ankle, leg and foot swelling.
- Reduce leg swelling during pregnancy.

**| Material |**

Lycra, Nylon.

Size 1 – 5 (Refer to P.29)

Color  



### Compression Pantyhose (Closed toe)

**| Application |**

- Minor varicosity.
- Tired, Aching Legs.
- Minor ankle, leg and foot swelling.

**| Material |**

Lycra, Nylon.

Size 1 – 5 (Refer to P.29)

Color  



Compression level  
**Class 2 (23-32 mmHg)**



**Knee high compression stockings (Closed toe)**

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



**Knee high compression stockings (Open toe)**

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



**Mid thigh compression stockings (Closed toe)**

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



**Mid thigh compression stockings (Open toe)**

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.
- Reduce leg swelling and prevent varicosity during pregnancy.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



# Compression stockings

Compression level

**Class 2 (23-32 mmHg)**



CPS-2302

## Thigh high compression stockings with silicone lace band (Closed toe)

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



CPS-2304

## Thigh high compression stockings with silicone lace band (Open toe)

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



CPS-2402

## Compression Pantyhose (Closed toe)


**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



CPS-2404

## Compression Pantyhose (Open toe)

**| Application |**

- Moderate to severe varicosity or spider veins.
- Post surgical use.
- Moderate edema.

**| Material |**

Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



Compression level  
**Class 3 (34-46 mmHg)**



**Knee high compression stockings (Closed toe)**

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

Size	1 – 7 (Refer to P.29)
Color	 

**| Material |**  
Lycra, Nylon.



**Knee high compression stockings (Open toe)**

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

Size	1 – 7 (Refer to P.29)
Color	 

**| Material |**  
Lycra, Nylon.



**Mid thigh compression stockings (Closed toe)**

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

Size	1 – 7 (Refer to P.29)
Color	 

**| Material |**  
Lycra, Nylon.



**Mid thigh compression stockings (Open toe)**

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

Size	1 – 7 (Refer to P.29)
Color	 

**| Material |**  
Lycra, Nylon.



# Compression stockings

Compression level

**Class 3 (34-46 mmHg)**



## Thigh high compression stockings with silicone lace band (Closed toe)

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

**| Material |**  
Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 




## Thigh high compression stockings with silicone lace band (Open toe)

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

**| Material |**  
Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 




## Compression Pantyhose (Closed toe)

- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

**| Material |**  
Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



## Compression Pantyhose (Open toe)


- | Application |**
- Severe varicosity.
  - Severe edema.
  - Post surgical use.

**| Material |**  
Lycra, Nylon.

Size	1 – 7 (Refer to P.29)
Color	 



**Size** | Compression stockings



Size		1	2	3	4	5	6	7
Circumference	cH	≤ 95 cm ≤ 37½"	≤ 100 cm ≤ 39¼"	≤ 105 cm ≤ 41¼"	≤ 110 cm ≤ 43¼"	≤ 114 cm ≤ 45"	≤ 118 cm ≤ 46½"	≤ 120 cm ≤ 47¼"
	cG	43 - 48 cm 17" - 18¾"	45 - 52 cm 17¾" - 19¾"	47 - 56 cm 18½" - 20½"	50 - 60 cm 19¾" - 22"	53 - 63 cm 20¾" - 22¾"	56 - 66 cm 22" - 23½"	58 - 68 cm 22¾" - 25¼"
	cC	28 - 34 cm 11" - 13¼"	30 - 37 cm 11¾" - 14½"	33 - 40 cm 13" - 15¾"	35 - 43 cm 13¾" - 17"	37 - 45 cm 14½" - 17¾"	39 - 47 cm 15¼" - 18½"	41 - 49 cm 16" - 19¼"
	cB	18 - 20 cm 7" - 7¾"	20 - 22 cm 7¾" - 8¾"	22 - 24 cm 8¾" - 9½"	24 - 26 cm 9½" - 10¼"	26 - 28 cm 10¼" - 11"	28 - 30 cm 11" - 11¾"	30 - 32 cm 11¾" - 12½"
Size		1	2	3	4	5	6	7
Length	IG	62 - 70 cm 24½" - 27½"	62 - 70 cm 24½" - 27½"	66 - 74 cm 26" - 29"	66 - 74 cm 26" - 29"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"
	ID	30 - 38 cm 11¾" - 15"	30 - 38 cm 11¾" - 15"	34 - 40 cm 13¼" - 15¾"	34 - 40 cm 13¼" - 15¾"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"

cH The circumference of hip. cG The circumference of thigh. cC The circumference of calf.  
cB The circumference of ankle. IG The length from heel to hip. ID The length from heel to knee.

## Diabetic socks



CPS-4001

### Diabetic socks

**| Feature |**


- Antibacterial silver fiber to prevent the growing of bacteria.
- Seamless stitching to prevent the friction on toes.
- Elastic and pressure free top band to avoid oppression.
- Comfortable toweling sole provides cushion and protection for feet.

**| Material |**

Combed cotton,  
Lycra covered yarn,  
Bamboo Charcoal fiber,  
Silver fiber.

**| Application |**

- Health caring for diabetic.
- People with bad circulation on foot.
- People who have foot odor.



SIZE	EURO	UK	US	
			MEN	WOMEN
M	38 - 42	4 - 7.5	4½ - 8	6 - 9½
L	43 - 47	6 - 10	6½ - 10½	8 - 11½

# Compression stockings

## Anti-embolism stocking



### Knee high anti-embolism stockings 18mmHg


#### | Application |

- For post-surgical, reclining or bed-confined patients to get better circulation.
- Reduce the risk of embolism formation.
- Prevent thrombo-embolic disease.

#### | Material |

Spandex, Nylon.

Size	1 - 7
Color	White



### Thigh high anti-embolism stockings with silicone top band 18mmHg

#### | Application |

- For post-surgical, reclining or bed-confined patients to get better circulation.
- Reduce the risk of embolism formation.
- Prevent thrombo-embolic disease.


#### | Material |

Spandex, Nylon.

Size	1 - 7
Color	White



### Size | Anti-embolism stockings



Size		1	2	3	4	5	6	7
Circumference	cH	≤ 95 cm ≤ 37½"	≤ 100 cm ≤ 39¼"	≤ 105 cm ≤ 41¼"	≤ 110 cm ≤ 43¼"	≤ 114 cm ≤ 45"	≤ 118 cm ≤ 46½"	≤ 120 cm ≤ 47¼"
	cG	43 - 48 cm 17" - 18¾"	45 - 50 cm 17¾" - 19¾"	47 - 52 cm 18½" - 20½"	50 - 56 cm 19¾" - 22"	53 - 58 cm 20¾" - 22¾"	56 - 60 cm 22" - 23½"	58 - 64 cm 22¾" - 25¼"
	cC	28 - 34 cm 11" - 13¼"	30 - 37 cm 11¾" - 14½"	33 - 40 cm 13" - 15¾"	35 - 43 cm 13¾" - 17"	37 - 45 cm 14½" - 17¾"	39 - 47 cm 15¼" - 18½"	41 - 49 cm 16" - 19¼"
	cB	18 - 20 cm 7" - 7¾"	20 - 22 cm 7¾" - 8¾"	22 - 24 cm 8¾" - 9½"	24 - 26 cm 9½" - 10¼"	26 - 28 cm 10¼" - 11"	28 - 30 cm 11" - 11¾"	30 - 32 cm 11¾" - 12½"
Length	IG	62 - 70 cm 24½" - 27½"	62 - 70 cm 24½" - 27½"	66 - 74 cm 26" - 29"	66 - 74 cm 26" - 29"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"	68 - 76 cm 26¾" - 30"
	ID	30 - 38 cm 11¾" - 15"	30 - 38 cm 11¾" - 15"	34 - 40 cm 13¼" - 15¾"	34 - 40 cm 13¼" - 15¾"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14¼" - 16½"

cH The circumference of hip. cG The circumference of thigh. cC The circumference of calf.  
cB The circumference of ankle. IG The length from heel to hip. ID The length from heel to knee.

## Travel stockings



CPS-1001

### Knee high travel stocking 10-14 mmHg

**| Application |**

- Minor ankle, leg and foot swelling during travel.
- Helps prevent DVT.

**| Material |**

Spandex, Nylon.

Size	One Size
Color	 



CPS-1002

### Knee high travel stocking 18-21 mmHg

**| Application |**

- Minor ankle, leg and foot swelling during travel.
- Helps prevent DVT.


**| Material |**

Spandex, Nylon.

Size	S - XL
Color	 



### Size | Travel stockings




Size		S	M	L	XL
Circumference	cC	28 - 34 cm 11" - 13¼"	33 - 40 cm 13" - 15¾"	37 - 45 cm 14½" - 17¾"	41 - 49 cm 16" - 19¼"
	cB	18 - 20 cm 7" - 7¾"	22 - 24 cm 8¾" - 9½"	26 - 28 cm 10¼" - 11"	30 - 32 cm 11¾" - 12½"
Size		S	M	L	XL
Length	ID	30 - 38 cm 11" - 15"	34 - 40 cm 13¼" - 15" <sup>W</sup>	36 - 42 cm 14¼" - 16½"	36 - 42 cm 14" <sup>W</sup> - 16" <sup>W</sup>

cC The circumference of calf. cB The circumference of ankle. ID The length from heel to knee.

# Insoles and Footcare series

Keep Comfort Health Plus



- Medical grade silicone barely cause allergies.
- Medical grade silicone never flattens, deforms or degrades for long time use.
- Gently massage the feet area.
- Help relieve foot pain and fatigue from prolonged standing or walking.
- Redistributes pressure on the foot.
- The ergonomic design of shape helps to easy on and off.
- Silicone and gel products are washable and reusable.



Medical grade silicone barely cause allergies.



Multiple density design helps to release and redistribute pressure on the foot.



The ergonomic design of shape helps to easy on and off.



Prevent uncomfortable friction.



Skin moisturizing.



Slim and light weight design.





# Insoles and Footcare

## Silicone



CPF-1001

### Silicone insoles

#### | Feature |

- Softer inserts absorb shock and pain on the metatarsal and heel areas.
- Raised metatarsal pad helps to support metatarsal arch.
- Multiple holes for air circulation.

SIZE	EURO	UK	US	
			MEN	WOMEN
34	33 - 34	1 - 2	2 - 3	3½ - 4½
36	35 - 36	2½ - 3½	3½ - 4½	5 - 6
38	37 - 38	4 - 5	5 - 6	6½ - 7½
40	39 - 40	5½ - 6½	6½ - 7½	8 - 9
42	41 - 42	7 - 8	8 - 9	9½ - 10
44	43 - 44	8½ - 9½	9½ - 10½	10½ - 11½
46	45 - 46	10 - 11	11 - 12	12 - 13½
48	47 - 48	11½ - 13	12½ - 13½	14 - 15



**| Material |**  
Silicone

**QTY : 1 PR**

#### | Application |

- Heel pain.
- Metatarsalgia.
- Relieves the discomfort of plantar fasciitis.




CPF-1002

### Silicone insoles with toe support

#### | Feature |

- Softer inserts absorb shock and pain on the lateral of mid foot and heel areas.
- Raised toe support helps to support toe area and relieve the discomfort of toes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	36 - 38	3½ - 5	4½ - 6	6 - 7½
M	39 - 41	5½ - 7½	6½ - 8½	8 - 9½
L	42 - 44	8 - 9½	9 - 10½	10 - 11½



**| Material |**  
Silicone

**QTY : 1 PR**

#### | Application |

- Heel pain.
- Relieves the pressure in lateral side.




CPF-1003

### Silicone insoles

#### | Feature |

- Soft liner inside to provide extra comfortable.
- Can be cut to fit every shoe size.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	34 - 37	2 - 4½	3 - 5½	4½ - 7
M	38 - 42	5 - 8	6 - 9	7½ - 10
L	43 - 47	8½ - 11½	9½ - 12½	10½ - 14



**| Material |**  
Silicone, Nylon

**QTY : 1 PR**

#### | Application |

- Foot pain.






CPF-1004

### Silicone insoles

**| Feature |**

- Softer inserts absorb shock and pain on the metatarsal and heel areas.
- Raised metatarsal pad helps to support metatarsal arch.

SIZE	EURO	UK	US	
			MEN	WOMEN
XS	37 - 38	4 - 5	5 - 6	6½ - 7½
S	39 - 40	5½ - 6½	6½ - 7½	8 - 9
M	41 - 42	7 - 8	8 - 9	9½ - 10
L	43 - 44	8½ - 9½	9½ - 10½	10½ - 11½
XL	45 - 46	10 - 11	11 - 12	12 - 13½



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel pain.
- Metatarsalgia.
- Relieves the discomfort of plantar fasciitis.



CPF-1007

### Silicone insoles

**| Feature |**

- Softer inserts absorb shock and pain on the arch and heel area.
- Slim and light weight design.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 40	4 - 6½	5 - 7½	6½ - 9
L	41 - 44	7 - 9½	8 - 10½	9½ - 11½



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel pain.
- Relieves the pressure in medial side.




CPF-1101

### Silicone insoles ¾ length

**| Feature |**

- Softer inserts absorb shock and pain on the metatarsal and heel areas.
- Provides heel stability and supports the longitudinal arch.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 40	4 - 6½	5 - 7½	6½ - 9
L	41 - 44	7 - 9½	8 - 10½	9½ - 11½



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel pain.
- Metatarsalgia.
- Relieves the discomfort of plantar fasciitis.



# Insoles and Footcare



CPF-1201

## Silicone heel cups

### | Feature |

- Softer insert provides more cushioning and shock absorption on heel area.
- Cupped design provides heel stability and comfort.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Heel Pain.
- Applicable for long time standing.




CPF-1202

## Silicone heel pads with removable pads

### | Feature |

- Helps absorb impact and prevent heel pain.
- Removable pads helps relieve heel spur.
- Could be use with/without removable pads for different needs.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Heel Pain.
- Heel spur.




CPF-1203

## Silicone heel pads

### | Feature |

- Softer insert provides more cushioning and shock absorption on heel area.
- Redistributes pressure on the foot.
- Non-slip design on the heel area.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Heel Pain.
- Applicable for long time standing.






CPF-1204

### Silicone heel cushions

**| Feature |**

- Softer insert on heel area helps absorb shock and relieve pain.
- Redistributes pressure on the foot.
- Different thickness of medial and lateral edge could be used for pronation or supination.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	≤ 34	≤ 2	≤ 3	≤ 4½
M	35 - 40	2½ - 6½	3½ - 7½	5 - 9
L	41 - 46	7 - 11	8 - 12	9½ - 13½



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel Pain.
- Applicable for long time standing.
- Slight varus or valgus.




CPF-1205

### Silicone heel pads

**| Feature |**

- Squares design helps absorb shock and relieve pain on heel area.
- Redistributes pressure on the foot.
- Thin type fits most shoe styles.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel Pain.
- Applicable for long time standing.



CPF-1206

### Silicone heel spur pads

**| Feature |**

- U-shaped design helps absorb shock and reduce pressure on the heel area.
- Provides heel stability and prevents ankle sprain.
- Relieves heel pain and heel spur.
- Easy to fit the most heel shape of the shoes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Heel Pain.
- Calcaneus instability.
- Flatfoot.



# Insoles and Footcare



CPF-1208

## Silicone heel cups

### | Feature |

- Soft material provides cushioning and shock absorption on heel area.
- Cupped design provides heel stability and prevents ankle sprain.
- Ergonomic design fits the shape of foot, provides comfortable wearing.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Heel Pain.
- Applicable for long time standing.



CPF-1301

## Silicone metatarsal pads

### | Feature |

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.
- Helps realign toes.
- Toe spreader design could be used for hallux valgus.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Metatarsalgia.
- Forefoot calluses.
- Overlapping toes.



CPF-1302

## Silicone metatarsal pads

### | Feature |

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Metatarsalgia.
- Forefoot calluses.





CPF-1401

### Silicone toe caps

**| Feature |**

- Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Fits for your toes or fingers.

S	Diameter 16mm
L	Diameter 20mm



**| Material |**  
Silicone

**QTY : 2 PCS**

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



CPF-1402

### Silicone toe tubes

**| Feature |**

- Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Separate and straighten overlapping toes.

S	Diameter 15mm
M	Diameter 18mm
L	Diameter 22mm



**| Material |**  
Silicone

**QTY : 2 PCS**

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



CPF-1501

### Silicone toe crests

**| Feature |**

- Helps reduce pressure and friction.
- Stretch contracted toes.
- Relieves pain caused by hammer, claw, and Mallet toes.
- Toe loop design keeps toe crest in place.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Hammer toes.
- Claw toes.
- Mallet toes.



# Insoles and Footcare



CPF-1502

## Silicone toe spreaders

### | Feature |

- Helps straighten and realign toes.
- Relieves pain caused by overlapping toes and hallux valgus.
- Reduces friction between toes.
- Separates the overlapping toes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½



| Material |  
Silicone

QTY : 2 PCS

### | Application |

- Hallux valgus.
- Overlapping toes.



CPF-1503

## Silicone toe separators

### | Feature |

- Helps straighten and realign toes.
- Relieves pain caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½



| Material |  
Silicone

QTY : 2 PCS

### | Application |

- Hallux valgus.
- Overlapping toes..



CPF-1504

## Silicone tailor's bunion protector

### | Feature |

- Helps reduce pressure and friction on 5th toe area.
- Toe loop design keeps protector in place.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



| Material |  
Silicone

QTY : 1 PC

### | Application |

- Tailor's bunion.
- Calluses.





### Silicone hallux valgus protector

**| Feature |**

- Helps reduce pressure and friction on hallux area.
- Toe loop design keeps protector in place.

**| Application |**

- Hallux valgus.
- Bunion.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



**| Material |**  
Silicone

**QTY : 1 PC**



### Silicone corn pad

**| Feature |**

- Soft silicone pad protects and prevents of corns and callous formation.

**| Application |**

- Calluses.
- Corns.

One Size



**| Material |**  
Silicone, Polyester

**QTY : 4 PCS**



### Silicone corn pads

**| Feature |**

- Soft silicone pad protects and prevents of corns and callous formation.

**| Application |**

- Calluses.
- Corns.

One Size



**| Material |**  
Silicone, Polyester

**QTY : 9 PCS**



### Silicone heel liners

**| Feature |**

- Prevents the friction between shoes and back of the heel.

**| Application |**

- Back of heel frictions.
- Calluses.

One Size



**| Material |**  
Silicone

**QTY : 1 PR**



# Insoles and Footcare

## Gel



CPF-2201

### Gel arch supports

**| Feature |**

- For arch support.
- Increases the stability of the foot.
- Prevents the arch pain.
- Fits most shoe styles.

One Size



**| Material |**  
Gel

**QTY : 1 PR**

**| Application |**

- Flatfoot.
- Over pronation foot.



CPF-2401

### Gel cushion

**| Feature |**

- Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Separate and straighten overlapping toes.
- Fully coated with gel.
- Could be cut according to the need.

S Width 22mm x 15cm

L Width 30mm x 15cm



**| Material |**  
Gel, Nylon, Spandex

**QTY : 1 TUBE**

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



CPF-2402

### Gel toe tubes

**| Feature |**

- Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Separates and straighten overlapping toes.

S Diameter 12mm

L Diameter 15mm



**| Material |**  
Gel

**QTY : 2 PCS**

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.





CPF-2403

### Gel toe caps

**| Feature |**

- Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Fits for your toes or fingers.
- Fully coated with gel.

S	Width 20mm
M	Width 25mm
L	Width 30mm



**| Material |**

Gel, Nylon, Spandex

**QTY : 2 PCS**

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



CPF-2404  
CPF-2405

### Gel toe pads CPF-2404

### Gel toe pads (roll) CPF-2405

**| Feature |**

- Gel dot provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Separates the overlapping toes.
- Could be cut according to the need. (CPF-2405)

	CPF-2404	CPF-2405
S	Width 20mm	Width 20mm x 55cm (17 dots)
M	Width 25mm	Width 25mm x 55cm (15 dots)
L	Width 30mm	Width 30mm x 55cm (11 dots)



**| Material |**

Gel, Nylon, Spandex

**QTY : 2 PCS**

(CPF-2404)

**QTY : 1 PC**

(CPF-2405)

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



CPF-2406  
CPF-2407

### Gel toe strip CPF-2406

### Gel toe strip (roll) CPF-2407

**| Feature |**

- Gel pad Provides cushioning and protection to toes.
- Relieves pain caused by calluses and corns.
- Separates the overlapping toes.
- Half coated with gel.
- Could be cut according to the need. (CPF-2407)

	CPF-2406	CPF-2407
S	Width 20mm	Width 20mm x 60cm
M	Width 25mm	Width 25mm x 60cm
L	Width 30mm	Width 30mm x 60cm



**| Material |**

Gel, Nylon, Spandex

**QTY : 2 PCS**

(CPF-2406)

**QTY : 1 PC**

(CPF-2407)

**| Application |**

- Overlapping toes.
- Calluses on toe.
- Corns on toe.



# Insoles and Footcare



CPF-2501

## Gel toe separators

### | Feature |

- Helps straighten and realign toes.
- Relieves pain caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



**| Material |**  
Gel

**QTY : L-2 PCS**  
**S-6 PCS**

### | Application |

- Hallux valgus.
- Overlapping toes.



CPF-2502

## Gel toe spreaders

### | Feature |

- Helps straighten and realign toes.
- Relieves pain caused by overlapping toes and hallux valgus.
- Reduces friction between toes.
- Separates the overlapping toes.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½



**| Material |**  
Gel

**QTY : 2 PCS**

### | Application |

- Hallux valgus.
- Overlapping toes.



CPF-2503

## Gel toe crests

### | Feature |

- Helps reduce pressure and friction.
- Stretch contracted toes.
- Relieves pain caused by hammer, claw, and Mallet toes.
- Toe loop design keeps toe crest in place.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 36	1 - 3½	2 - 4½	3½ - 6
M	37 - 41	4 - 7	5 - 8	6½ - 9½
L	42 - 46	7½ - 11	8½ - 12	10 - 13½



**| Material |**  
Gel

**QTY : 1 PR**

### | Application |

- Hammer toes.
- Claw toes.
- Mallet toes.





### Gel bunion protector

**| Feature |**

- Helps reduce pressure and friction on hallux area.
- Toe loop design keeps protector in place.

**| Application |**

- Hallux valgus.
- Bunion.

One Size



**| Material |**  
Gel

**QTY : 1 PC**



### Gel bunion protector with toe spreader

**| Feature |**

- Helps straighten and realign toes.
- Relieves pain caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

**| Application |**

- Hallux valgus.
- Bunion.
- Overlapping toes.

One Size



**| Material |**  
Gel

**QTY : 1 PC**



### Gel toe spreaders

**| Feature |**

- Relieves pain caused by overlapping toes and corns.
- Reduces friction between toes.
- Separates the overlapping toes.

**| Application |**

- Hallux valgus.
- Bunion.
- Overlapping toes.

One Size



**| Material |**  
Gel, Nylon, Spandex

**QTY : 1 PR**



### Metatarsal foot pads

**| Feature |**

- Gel pad provides more cushioning and shock absorption on metatarsal area.
- Elastic sleeve design keeps metatarsal foot pad in place.

**| Application |**

- Metatarsalgia.
- Forefoot calluses.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	33 - 38	1 - 5	2 - 6	3½ - 7½
L	39 - 46	5½ - 11	6½ - 12	8 - 13½



**| Material |**  
Gel, Nylon, Spandex

**QTY : 1 PR**



# Insoles and Footcare



CPF-6302

## Metatarsal wrap

### | Feature |

- Gel pad provides more cushioning and shock absorption on ball of foot area.
- Prevents calluses and relieve pain on metatarsal area.
- Velcro design for easy adjusting.

### | Application |

- Metatarsalgia.
- Forefoot calluses.

One Size



| Material |  
Gel, PVC, Polyester

QTY : 1 PC



CPF-6601

## Gel bunion sleeve

### | Feature |

- Gel pad provides more cushioning and shock absorption on hallux area.
- Helps realign toes.

### | Application |

- Hallux valgus.
- Bunion.

SIZE	EURO	UK	US	
			MEN	WOMEN
S	35 - 37	2½ - 4	3½ - 5	5 - 6½
M	38 - 40	4½ - 6½	5½ - 7½	7 - 9
L	41 - 43	7 - 8½	8 - 9½	9½ - 11



| Material |  
Gel, Nylon, Spandex

QTY : 1 PC



CPF-6602

## Achilles protection pad

### | Feature |

- Relieves pressure and pain caused by Achilles tendonitis.
- Soft gel pad helps reduce friction and absorb shock on the Achilles area.

### | Application |

- Achilles tendonitis.

S	13 - 17 cm (5" - 6¾")
L	17 - 25 cm (6¾" - 9¾")



| Material |  
Gel, Nylon, Spandex

QTY : 1 PC



CPF-6604

## Gel heel socks

### | Feature |

- Relieves pressure and pain around heel area.
- Gel coated on heel area.
- Mineral oil helps soften and moisturize skin.

### | Application |

- Calluses and stratum corneum.
- Dry skin.

One Size



| Material |  
Gel, Nylon, Spandex

QTY : 1 PR



## High heel & Sandals



### High heel silicone insoles

**| Feature |**

- Raised metatarsal pad helps to support metatarsal arch.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

SIZE	EURO	UK	US
S	35 - 37	2½ - 4	5 - 6½
M	38 - 40	4½ - 6½	7 - 9
L	41 - 43	7 - 8½	9½ - 11



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Suitable for high heel shoes.
- Metatarsalgia.



### High heel silicone insoles

**| Feature |**

- Softer inserts absorb shock and pain on the metatarsal area.
- Raised metatarsal pad helps to support metatarsal arch.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

SIZE	EURO	UK	US
S	35 - 37	2½ - 4	5 - 6½
M	38 - 40	4½ - 6½	7 - 9
L	41 - 43	7 - 8½	9½ - 11



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Suitable for high heel shoes.
- Metatarsalgia.



### High heel Silicone heel pads

**| Feature |**

- Helps absorb shock and relieve pain on heel area.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

SIZE	EURO	UK	US
S	33 - 37	1 - 4	3½ - 6½
L	38 - 42	5 - 8	7 - 10



**| Material |**  
Silicone

**QTY : 1 PR**

**| Application |**

- Suitable for high heel shoes.
- Heel pain.



# Insoles and Footcare



CPF-1303

## High heel Silicone metatarsal pads

### | Feature |

- Raised metatarsal pad helps to support metatarsal arch.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

SIZE	EURO	UK	US
S	33 – 37	1 – 4	3½ - 6½
L	38 – 42	5 – 8	7 – 10



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Suitable for high heel shoes.
- Metatarsalgia.



CPF-1304

## High heel Silicone metatarsal pads with removable dot

### | Feature |

- Raised metatarsal pad helps to support metatarsal arch.
- Removable dot could be adjusted for different support region.
- Slim and light weight design fits most women's fashion footwear and high heel shoes style.

SIZE	EURO	UK	US
S	33 – 37	1 – 4	3½ - 6½
L	38 – 42	5 – 8	7 – 10



| Material |  
Silicone

QTY : 1 PR

### | Application |

- Suitable for high heel shoes.
- Metatarsalgia.



CPF-2301

## Gel metatarsal pads

### | Feature |

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.
- Soft textile makes comfortable hand feel.
- Thin design for most high heel and casual shoes.

One Size



| Material |  
Gel, Polyester

QTY : 1 PR

### | Application |

- Suitable for high heel shoes.
- Metatarsalgia.
- Forefoot calluses.





### Gel metatarsal pads

**| Feature |**

- Provides more cushioning and shock absorption on metatarsal area.
- Reduces and prevents pain caused by calluses.
- Thin design for most high heel and casual shoes.

One Size



**| Material |**  
Gel

**QTY : 1 PR**

**| Application |**

- Suitable for high heel shoes.
- Metatarsalgia.
- Forefoot calluses.



### Metatarsal pads for sandals

**| Feature |**

- Prevents the friction between toes and sandals.
- Reduces the pain caused by calluses on the metatarsal area.
- Relieve the pressure and friction on metatarsal area during long-time standing and walking.
- The adhesive bottom is convenient to plaster on the surface of the sandals and it is immovable.

One Size



**| Material |**  
PU gel

**QTY : 1 PR**

**| Application |**

- Suitable for sandals.
- Metatarsalgia.
- Forefoot calluses.
- Toe friction.



### Toe protector for sandals

**| Feature |**

- Prevents the friction between toes and sandals.
- Design for pinching protection, and can be use direct on flip sandals.
- The adhesive bottom is convenient to plaster on the surface of the sandals and it is immovable.

One Size



**| Material |**  
PU gel

**QTY : 1 PR**

**| Application |**

- Suitable for sandals.
- Toe friction.



# Insoles and Footcare

## Moisturizing



### Moisturizing socks

#### | Feature |

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
- Moisturizing gel improves dry skin.
- Lavender essential oil helps recovering weary spirit.

One Size

#### | Material |

Outer sleeve: Nylon, Spandex.  
Inner gel: Gel, Jojoba oil,  
Lavender essential oil.

QTY : 1 PR

#### | Application |

- Calluses and stratum corneum.
- Dry skin.



### Moisturizing heel socks

#### | Feature |

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
- Moisturizing gel improves dry skin.
- Lavender essential oil helps recovering weary spirit.

One Size

#### | Material |

Outer sleeve: Nylon, Spandex.  
Inner gel: Gel, Jojoba oil,  
Lavender essential oil.

QTY : 1 PR

#### | Application |

- Calluses and stratum corneum.
- Dry skin.



### Moisturizing gloves

#### | Feature |

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
- Moisturizing gel improves dry skin.
- Lavender essential oil helps recovering weary spirit.

One Size

#### | Material |

Outer sleeve: Nylon, Spandex.  
Inner gel: Gel, Jojoba oil,  
Lavender essential oil.

QTY : 1 PC

#### | Application |

- Calluses and stratum corneum.
- Dry skin.



### Moisturizing elbow sleeves

#### | Feature |

- Sustained release jojoba oil helps to soften calluses and stratum corneum.
- Moisturizing gel improves dry skin.
- Lavender essential oil helps recovering weary spirit.

One Size

#### | Material |

Outer sleeve: Nylon, Spandex.  
Inner gel: Gel, Jojoba oil,  
Lavender essential oil.

QTY : 1 PC

#### | Application |

- Calluses and stratum corneum.
- Dry skin.



**COMFORT PLUS CORPORATION**

ADD | 14F-1, No.96, Xingde Rd., Sanchong Dist., New Taipei City 24158, Taiwan

TEL: +886-2-85113758 FAX: +886-2-85113755 E-mail | [service@comfortplus-medical.com](mailto:service@comfortplus-medical.com)

<http://www.comfortplus-medical.com>